

Door Handles and Adapters for Omnicell® FlexLock

Special accessories are available for facilities that install the Omnicell FlexLock on Helmer upright refrigerators, undercounter refrigerators, and undercounter freezers.

FlexLock Compatible Door Handles (Upright Models)

The door handle features a cutout and spacer which enable the FlexLock to wrap around the front and side of the door.

- Available for single door upright refrigerators 11, 20, 25, and 26.5 cu ft (326, 572, 714 and 750 l)
- Compatible with glass and solid doors
- Factory or field installation


FlexLock Adapter Kit (Undercounter Models)

The adapter kit includes an extension bar that mounts to the door of the unit, allowing the FlexLock to be properly installed.

- Adapter kits available for undercounter refrigerators and freezers (5 cu ft / 142 l)
- Compatible with glass and solid doors; powder coated and stainless steel exteriors
- Factory installed or field installation kits available



Cut-out of Door Handle For Upright Refrigerators



Upright Omnicell®
FlexLock Compatible
Door Handle



Undercounter Adapter
(Shown with Omnicell®
FlexLock)

Part No.	Description
400848-2	Omnicell FlexLock Adapter Kit , Factory installed, 5 cf (142 l) (for UC Refrigerators and Freezers)
400782-2	Omnicell FlexLock Compatible Door Handle , Factory installed, 11 cf (326 l)
400782-1	Omnicell FlexLock Compatible Door Handle , Factory installed, 20, 25 and 26.5 cf (572, 714 and 750 l)
400848-1	Omnicell FlexLock Adapter Kit , Field installation, 5 cf (142 l) (for UC Refrigerators and Freezers)
400875-2	Omnicell FlexLock Compatible Door Handle , Field installation, 11 cf (326 l)
400875-1	Omnicell FlexLock Compatible Door Handle , Field installation, 20, 25 and 26.5 cf (572, 714 and 750 l)

Note: Adapter kits and compatible door handles do not include the remote lock