


Accessories


Bridge Shelf for Refrigerators

The Bridge Shelf spans the space between two shelves in a Helmer double door refrigerator offering over 13% more storage area per shelf.

- Available for 45 cu ft and 56 cu ft models
- Attaches to existing shelves (no additional hardware required)
- Bacteria-resistant powder coated



45 cu ft Bridge Shelf (22" w x 18" d)



56 cu ft Bridge Shelf (22" w x 24" d)

Part # 400845-1 (18")
Bridge shelf for 45 cu ft refrigerator

400845-2 (24")
Bridge shelf for 56 cu ft refrigerator

Increased Shelf Space

| | 45 cu ft models | 56 cu ft models |
|----------------------|----------------------|----------------------|
| Without Bridge Shelf | 968 in ² | 1276 in ² |
| With Bridge Shelf | 1100 in ² | 1452 in ² |
| Additional Space | 13.6% | 13.8% |



www.helmerinc.com
800.743.5637